

Instructions For Use

Always wash your hands thoroughly before connecting any Extended Wear Bag.

Maintenance with Catheters

- The Extended Wear Leg Bag is designed to maintain a "closed" system with catheter use
- It is very important that a sterile technique is required for connection to catheters
- A new sterile leg bag and tubing should be attached when a new catheter is inserted
- Do not disconnect any bag from the catheter until the catheter is due for a change or this should be done at the discretion of your Health Professional.

Connecting the Leg Bag

1. Cut a small length of the extension tubing and attach to the leg bag outlet to facilitate night drainage bag connection.
2. Cut remaining tubing to desired length, connect one end to the leg bag inlet and the other end to the sheath or catheter.
3. Attach leg bag straps to accommodate horizontal or vertical placement on thigh or calf.
4. Open tap to 90 degrees for emptying.
5. Loosen or undo leg bag straps for sleeping.
6. Tap may be turned to the side when open and attached to night drainage.

Sheath Drainage with Leg bag & connecting Night Drainage Bag

1. Cut a short length of the extension tubing and place on Leg Bag outlet.
2. Connect the remaining tubing to the Leg Bag inlet.
3. Connect extension tubing to penile sheath.
4. Attach Leg Bag Straps to accommodate either horizontal or vertical placement on thigh or calf.

For more information or samples please call Customer Engagement on
Australia Freecall: 1800 219 179
New Zealand Freecall: 0800 678 669

5. For night drainage attach the Night Drainage Bag to end of Leg Bag outlet.
6. Loosen or undo the Leg Bag Straps for sleeping, or wear in the horizontal position.
7. Remove Leg Bag and extension tubing daily to clean and dry.
8. Replace Leg Bag at sign of wear and tear. Extended wear times can be up to 3 months.

Cleaning the Extended Wear Leg Bag (for sheath drainage)

1. Change leg bag daily.
2. Empty all urine from leg bag into toilet.
3. Close outlet tap, fill leg bag with cold tap water and shake vigorously.
4. Empty water through drainage outlet into toilet and repeat steps 3 and 4 from Connecting the leg bag.
5. Disconnect tubing and pour 30mls decontaminant solution into ends of tubing and all leg bag outlets.
6. Reconnect tubing and close leg bag outlet. Pour remaining solution into leg bag and agitate for 30 seconds.
7. Drain solution into toilet, open outlet and hang leg bag in shaded place to air dry. All parts must be completely dry before the next use.
8. Replace leg bag at signs of wear and tear. Extended wear time of up to 3 months is common.

Decontaminant Solution Preparation (for sheath drainage)

- Option 1:** Household detergent solution
- Dilute 5ml detergent in 150ml of warm tap water
- Option 2:** Weak vinegar solution (for odour control)
- Dilute 30ml vinegar in 120ml of cold tap water
- Option 3:** Milton's Solution

Hollister Continence Care

Extended Wear Bags

People First.



Attention to Detail. Attention to Life.




Extended Wear Bags & Combination Packs

- Anti reflux valve
- Low profile buttons on the leg bag easily secure the wide, machine washable straps
- Bag can be worn vertically or horizontally for increased discretion
- Easy to open 90° lever tap is designed for people with limited or poor dexterity
- Universal tapered adaptor provides an easy connection to ensure maximum security
- Available as combination packs or as individual items

Extended Wear Pack includes:

- 1 x Leg bag
- 1 x Latex-free extension tubing 42cm
- 1 x Pair fabric leg straps -medium 58cm long

*Leg bag, straps and tubing are also sold separately



Made from durable polyvinyl, the Hollister Extended Wear Range offers client convenience and cost savings whilst reducing the risk of introducing an infection by maintaining a closed system.


Extended Wear: Combination Packs & Leg Bags		
Stock Number	Bag Size/Length	QTY
9348 Combination Pack	550ml	Box of 12
9349 Combination Pack	900ml	Box of 12
9814/10 Bag only	550ml	Box of 10
9805/10 Bag only	900ml	Box of 10

Extended Wear: Leg Straps		
9342/3	38cm	1 Pair
9343/3	58cm	1 Pair

Non Latex Extension tubing		
9346	42cm	Box of 10

Extended Wear Night Drainage Bag

- 2000ml capacity
- Wide bore drain tube clamp outlet, with tuck away pocket
- 120cm soft non-kink, wide bore inlet tubing- allows freedom of movement
- Universal tapered adaptor provides an easy connection to ensure maximum security
- Connected bed hanger and sheet clip-for added security
- Protected and self sealing sample ports
- Sterile fluid pathway
- Graduations-enable easy visual assessment of urine output



Extended Wear: Night Drainage Bag		
Stock Number	Bag Size	QTY
9839/20	2000ml	Box of 20