

# Faecal Collector

## Guide for Use

- Easy to apply odour-barrier pouch
- Flexend Skin barrier helps protect perianal skin from irritation and maceration
- Ultra thin barrier conforms to patients anatomy, moving and stretching with the patient.
- Fabric backing on the barrier improves flexibility and barrier cohesion
- Attaches to collection devices for continuous drainage
- Clear film
- Not made with natural rubber latex
- Wear time of up to 7 days



### Faecal Collector

Stock No.	Description	Box QTY
9823/10	500 ml	Box of 10


Application instructions on reverse.

# Faecal Collector

## Guide for Use

### Application

- Turn patient to one side, with upper knee slightly flexed
- Clean and dry perianal skin, remove excess hair if necessary
- Cut the opening in the skin barrier to fit patient's anatomy
- For female patients, trim inside the shaded perianal triangle area
- Remove the release paper and fold barrier in half
- Separate the buttocks and position skin barrier over the patient's anus
- Apply the barrier, hold in place for 30 seconds to activate adhesive bond
- Connect to a bedside drainage collector, or close the drain cap



### Removal

- Disconnect the drainage collector and close the cap
- Using the fabric tab, gently peel the skin barrier from the skin
- Fold the barrier in half and discard

### Application tips

- Remove excess hair for maximum adhesion
- Ensure patient's skin is completely dry and free from powders, lotions, ointments or oily residue
- Do not use skin gel or other skin sealants as they may decrease wear time
- To enhance the seal, Hollister Adapt Paste may be used around the opening in the skin barrier
- Adapt CeraRing Seals can also be used to smooth skin surfaces, to prevent faecal matter from channeling under the skin barrier
- Hollister Adapt Powder may be used on eroded perianal skin. Excess powder must be brushed away before applying the skin barrier.

Prior to use, be sure to read the Instructions for Use for information regarding Intended Use, Contraindications, Warnings, Precautions, and Instructions.

**For more information or samples please call:**

**Australia Customer Care 1800 880 851**  
**www.hollister.com.au**

**New Zealand Customer Service 0800 678 669**  
**www.hollister.co.nz**



To view a video of "Instructions For Use", please visit our website

Australia: [www.hollister.com.au](http://www.hollister.com.au) OR  
New Zealand: [www.hollister.co.nz](http://www.hollister.co.nz)

Search with keyword "faecal collector" and click on link to product page, scroll down for video instructions