


Proven to prevent the transmission of pathogens that may cause UTIs


The American Society for Testing and Materials (ASTM) is an international standards organisation that develops and publishes voluntary consensus technical standards for a wide range of materials. ASTM has also conducted and published approximately 6000 test methods to date.

ASTM F1671 is a pass or fail test designed to determine if protective materials prevent the transmission of blood-borne pathogens.


VaPro catheter sleeve material was tested. In the test, one side of the sleeve material was exposed to a solution containing a type of virus called bacteriophage. Results showed that the viruses did not pass through the sleeve material.




Viruses
(smaller in size)



Bacteria



Fungi



Parasites
(greater in size)

The VaPro catheter sleeve material passes the ASTM F1671 test, providing assurance that it protects against pathogens that may cause UTIs.

Go to Hollister.com.au or Hollister.co.nz to learn more about clinical and laboratory evidence that demonstrates the benefits of VaPro catheters' protective tip and sleeve.

Hollister Nurse Portal

Better Resources, Better Outcomes

Time saving solutions lead to improved patient outcomes

A revolution in sample ordering



Order combinations of Samples, Continence Care Kits and tailor the content and quantity to your liking.



Search by product specifications or code to quickly find your desired product, and access past orders for quick re-ordering of favourites.



The portal is open 24/7 for quick ordering from your desktop or smart phone at your convenience. Instantly download literature for top-quality patient education at the click of a button.



You will receive an email confirmation and record of your order, and our Customer Care team will process it within 1 business day.

Log in or Register for complete access to our Continence Nurse Portal
www.hollisternurseportal.com.au



Scan and sign in today!

Hollister Australia
PO Box 375
Box Hill, VIC 3128
www.hollister.com.au

Nurse Priority Line
1800 880 840

Hollister New Zealand
PO Box 62-027
Mt Wellington, Auckland
www.hollister.co.nz

Customer Service
0800 678 669



With 100% No Touch Protection,
whatever you touch
won't touch VaPro catheters.

VaPro catheters provide 100% No Touch Protection with a protective tip and sleeve

VaPro catheters are hygienic by design to help protect patients from germs throughout the entire catheterisation process.

Protective tip

- helps shield the sterile catheter during insertion from bacteria located within the first 15 mm of the distal urethra
- helps reduce the risk of carrying bacteria into the urinary tract




Protective sleeve


- allows for catheter to be gripped anywhere
- provides a barrier to the catheter to help keep germs away

The protective tip is designed for male and female use

The protective tip works to greatly reduce the risk of bacteria being carried into the urinary tract during the insertion process. This illustration demonstrates its effectiveness using the male anatomy.



The catheter is advanced into the protective tip, and then the protective tip is inserted into the urethra.



The protective tip helps prevent the catheter from coming into contact with the bacteria that are predominantly in the first 15 mm of the urethra.



The protective tip helps protect the catheter from becoming contaminated.

Watch the Contamination and Prevention in IC Treatment Video and learn more about Hollister tip & sleeve protection and how they are designed to help protect you from germs. Scan the QR code or visit bit.ly/NoTouchAU or bit.ly/NoTouchNZ



Australia



New Zealand

A recent study confirms high patient satisfaction

Study explored real-world preferences associated with using VaPro catheters

25 participants from 5 countries (UK, Germany, France, the Netherlands, and the US) were interviewed to assess the perceived advantages of VaPro catheters compared to other intermittent catheters.

All interviewees had experiences with competitive intermittent catheters and changed to VaPro catheters.



Participants reported a high level of **satisfaction** with VaPro catheters enabling them to overcome issues they had encountered with other intermittent catheters.



Some participants reported that VaPro catheters had a **positive** effect on their quality of life.



Participants find VaPro catheters to be **hygienic** and **easy** to use.




50% of participants reported **less** frequent UTIs after changing to VaPro catheters.

Easy to use

The combination of protective tip and sleeve creates a no touch system that is easy to use
Ready to use right out of the package with no extra steps required



Designed with smooth rounded eyelets to enhance user comfort during insertion and withdrawal



Packaging designed with finger holes to be easy to open



Pocket-sized packaging for easy transportation and discreet out of home usage



Australia
www.hollister.com.au/vapro



New Zealand
www.hollister.co.nz/vapro

To learn more about what users are saying about VaPro catheters, visit our website.

VaPro catheters: 100% No Touch Protection



VaPro™ No Touch Intermittent Catheters				
Stock Code	Colour Code	Size	Length	System
72104	● Black	10 Ch	40 cm	Nelaton
72124	○ White	12 Ch	40 cm	Nelaton
72144	● Green	14 Ch	40 cm	Nelaton
72164	● Orange	16 Ch	40 cm	Nelaton
72102	● Black	10 Ch	20 cm	Nelaton
72122	○ White	12 Ch	20 cm	Nelaton
72142	● Green	14 Ch	20 cm	Nelaton



VaPro Pocket™ No Touch Intermittent Catheter				
Stock Code	Colour Code	Size	Length	System
70104	● Black	10 Fr	40 cm	Nelaton
70124	○ White	12 Fr	40 cm	Nelaton
70144	● Green	14 Fr	40 cm	Nelaton
70164	● Orange	16 Fr	40 cm	Nelaton



VaPro Plus Pocket™ No Touch Intermittent Catheter with integrated Collection Bag				
Stock Code	Colour Code	Size	Length	System
71084	● Blue	8 Fr	40 cm	Nelaton
71104	● Black	10 Fr	40 cm	Nelaton
71124	○ White	12 Fr	40 cm	Nelaton
71144	● Green	14 Fr	40 cm	Nelaton
71164	● Orange	16 Fr	40 cm	Nelaton



Request a sample



Australia
www.hollister.com.au/vapro



New Zealand
www.hollister.co.nz/vapro

VaPro™
No Touch Intermittent Catheter