





Product	Specs/SAS Qtys	Description
Medical Adhesive Remover Spray 	7737 76 g Spray Can Each SAS Monthly issue 2	360 degree 'Bag-in-a-can' spray technology - can be sprayed at any angle with minimal product wastage. Eases removal of skin barriers or adhesive borders (tapes), from the skin. <ul style="list-style-type: none"> • Breaks down adhesive bonds for easy removal of barrier adhesives including tapes. • Helps reduce the potential for skin stripping.
Medical Adhesive Spray 	7730 90 g Spray Can SAS Monthly issue 1	A silicone adhesive available in an aerosol can. Used to increase the adhesion of skin barriers and/or adhesive border products.* <ul style="list-style-type: none"> • Helps improve the adhesion of skin barriers. • Helps to provide a secure seal. *Requires STN authorisation
Skin Cleanser 	7210 Pump Spray Bottle 236 ml SAS Monthly issue 1	<ul style="list-style-type: none"> • Gentle foaming cleanser in a convenient pump spray. • Helps maintain healthy skin condition.
Skin Conditioning Creme 	7220 Conditioning Creme Bottle 118 ml SAS Monthly issue 2	<ul style="list-style-type: none"> • Smooths cracked and roughened skin • Non greasy, non perfumed. • Conditions and soothes skin. • Designed for use with ostomy products. • Apply minimally.



Hollister Ostomy Accessories

The options you need.
The expertise you want.

For more information or samples please call Customer Engagement on:

Australia **1800 219 179**
New Zealand **0800 678 669**


www.hollister.com





Hollister Ostomy. **Details Matter.**[™]





Product	Specs/SAS Qtys	Description
---------	----------------	-------------

	7805 Small 48 mm Standard thickness	Flat, circular-shaped hydrocolloid skin barriers with a centre opening. Fill uneven skin surfaces around the stoma. Easily shaped to the dimensions of the stoma. <ul style="list-style-type: none"> • Improve skin barrier fit. • Minimise leakage risk. • Protect the skin. • Improve pouch wear time. • Can be stretched, moulded or stacked to improve the fit of the skin barrier. • Individually packaged.
	7815 Small 48mm Slim	
	7806 Large 98mm Slim Box of 10 SAS Monthly issue 30	


	79520 20 mm	A ring of hydrocolloid material with a flat side and an outward-facing curved (convex) side, with a circular opening. Creates a seal that can be stretched to the dimensions of the stoma. Provides gentle convexity to the skin around the stoma. <ul style="list-style-type: none"> • Improve skin barrier fit. • Minimise leakage risk. • Helps allow output from a flush or recessed stoma to drain into the pouch.
	79530 30 mm	
	79540 40 mm Box of 10 SAS Monthly issue 30	

	79601 22 x 38 mm	A ring of hydrocolloid material with a flat side and an outward-facing curved (convex) side, with an oval opening. Creates a seal that can be stretched to the dimensions of the stoma. Provides gentle convexity to the skin around the stoma. <ul style="list-style-type: none"> • Improve skin barrier fit. • Minimise leakage risk. • Helps allow output from a flush or recessed stoma to drain into the pouch.
	79602 30 x 48 mm	
	79603 38 x 56 mm Box of 10 SAS Monthly issue 30	


	79400 Box of 10 SAS Monthly issue 30	A blend of hydrocolloids in strip format. Can be stretched and moulded to create custom shapes. Fill that narrow cavity or crease with ease.
--	---	--


	7906 28 gm puff bottle SAS Monthly issue 2	A blend of different hydrocolloids in a powder mixture. Absorbs excess moisture from skin around the stoma. <ul style="list-style-type: none"> • Helps protect and heal raw and weeping skin. • Helps the skin barrier adhere more effectively to most skin.
--	---	--


Product	Specs/SAS Qtys	Description
---------	----------------	-------------

	79300 57gm tube SAS Monthly issue 3	A blend of hydrocolloids in a paste formulation. Fill any uneven surface. Low alcohol paste helps minimise stinging. <ul style="list-style-type: none"> • Improve skin barrier fit. • Minimise leakage risk. • Protects skin
	79301 14.2 gm tube Not on the SAS	

	78500 236 ml bottle SAS Monthly issue 1 bottle	<ul style="list-style-type: none"> • A lubricant and deodorant in the one product. • Neutralises odour and lubricates the inside surface of the pouch to limit pancaking and ease pouch emptying
	78501 8 ml sachet SAS Monthly issue 100 sachets per year (additional to bottle issue)	

	7300 66 - 109 cm Medium	<ul style="list-style-type: none"> • An elasticised, beige, fabric belt specifically designed to attach to an ostomy pouch. • Attaches to the pouch and wraps around the body. • Helps to keep the pouch snugly in place.
	7299 74 - 125 cm Large	
	Each Annual issue 4	

	7917 50 sachets per box SAS Monthly issue 2 boxes	<ul style="list-style-type: none"> • Helps prevent skin stripping resulting from the frequent removal of adhesive products or skin barriers. • Helps protect skin from ostomy output.
--	--	---

	7760 50 sachets per box SAS Monthly issue 2 boxes	Gentle removal of adhesive products and adhesive residue. <ul style="list-style-type: none"> • Makes removal of pouching systems fast and easy. • Helps reduce the potential for skin stripping.
---	--	--