

VaPro Plus Pocket
 Touch Free Hydrophilic Intermittent Catheter

Hollister Code	Length (cm)	Size	Box Qty
71104	40	10Fr	30
71124	40	12Fr	30
71144	40	14Fr	30
71164	40	16Fr	30



For more information or to request a product sample, visit hollister.com.au/vapropluspocket



For more information regarding Continence Care, visit our eLearning Centre hollister.com.au/continencecarelearningcentre

Hollister Customer Care
 1800 880 851
www.hollister.com.au

Prior to use of VaPro catheters, be sure to read the Instructions for Use for information regarding Intended Use, Contraindications, Warnings, Precautions, and Instructions.

The Hollister logo, VaPro and VaPro Plus Pocket are trademarks of Hollister Incorporated. All other trademarks and copyrights are the property of their respective owners. Not all products are CE marked. © 2020 Hollister Incorporated. AUHC233. January 2020.



With 100% No Touch Protection, whatever you touch won't touch VaPro Plus Pocket catheters.



Continence Care



Continence Care



VaPro Plus Pocket is offered in pocket-sized packaging for easy transportation and discreet out of home usage

VaPro Plus Pocket catheters provide 100% No Touch Protection with a protective tip and sleeve

No-touch catheter with an integrated collection bag in pocket-sized packaging

Protective Tip

- Helps shield the sterile catheter during insertion from bacteria located in the first 15 mm of the distal urethra.
- Help reduce the risk of carrying bacteria into the urinary tract.

Protective Sleeve

- Allows for catheter to be gripped anywhere
- Provide a barrier to the catheter to help keep germs away

Collection Bag

- No toilet or receptacle is needed with the integrated collection bag of VaPro Plus Pocket catheters
- The collection bag includes a handle for convenient handling when emptying the collection bag and anti-reflux valve to prevent backflow of urine.

