

## Advance Plus & Advance Plus Pocket: Pre-Lubricated Intermittent Catheter Closed Systems

Advance Plus pre-lubricated intermittent catheters have all the same features as Advance, with the added benefit of an integrated 1500 ml collection bag.

### Advance Plus Pocket

All the features of Advance Plus and...

- 60% smaller packaging.
- Easy to store and carry, making your daily routine more convenient.



ADVANCE  
PLUS POCKET:  
CONVENIENT  
& DISCREET

### Advance Plus: Pre-Lubricated Intermittent Catheters

Stock Code	Funnel Colour	Size/Gender	Length	Description	Box QTY
94104	●	10Fr Male	40 cm	Straight	25
94124	○	12Fr Male	40 cm	Straight	25
94144	●	14Fr Male	40 cm	Straight	25
94164	●	16Fr Male	40 cm	Straight	25
94184	●	18Fr Male	40 cm	Straight	25
94102	●	10Fr Female	20 cm	Straight	25
94122	○	12Fr Female	20 cm	Straight	25
94142	●	14Fr Female	20 cm	Straight	25

### Advance Plus Pocket: Pre-Lubricated Intermittent Catheters

Stock Code	Funnel Colour	Size	Length	Description	Box QTY
93104	●	10Fr	40 cm	Straight Pocket	25
93124	○	12Fr	40 cm	Straight Pocket	25
93144	●	14Fr	40 cm	Straight Pocket	25
93164	●	16Fr	40 cm	Straight Pocket	25

### Advance Plus Tiemann Tip: Pre-Lubricated Intermittent Catheters

Stock Code	Funnel Colour	Size	Length	Description	Box QTY
95124	○	12Fr	40 cm	Tiemann Tip	25
95144	●	14Fr	40 cm	Tiemann Tip	25
95164	●	16Fr	40 cm	Tiemann Tip	25

#### Australia

6/345 Ingles Street, Port Melbourne,  
Victoria, Australia 3207

#### New Zealand

58 Richard Pearse Drive, Airport Oaks,  
Auckland, New Zealand

[www.hollister.com](http://www.hollister.com)

Customer Engagement

Australia **1800 219 179**

New Zealand **0800 678 669**



HL033

## Hollister Continence Care



Intermittent Catheter Range

## Apogee: Intermittent Catheters

Apogee intermittent catheters are latex-free and ready to use with the addition of lubrication.

- **Smooth eyelets** provide comfort on insertion and removal
- Smooth 100% **medical grade PVC**
- **Easy to open** packaging
- **Full range** of sizes from 6Fr to 18Fr
- **Firm but flexible** for greater comfort and control

Available in both curved and straight packaging for convenience.




### Apogee: Intermittent Catheters

Stock Code	Funnel Colour	Size/Gender	Length	Description	Box QTY
1007	●	6Fr Male	38 cm	Firm	30
1018	●	8Fr Male	40 cm	Firm	30
1034	●	10Fr Male	40 cm	Firm	30
1046	○	12Fr Male	40 cm	<b>Soft</b>	30
1048	○	12Fr Male	40 cm	Firm	30
1049	○	12Fr Male	40 cm	Firm - Curved Packaging	<b>50</b>
1061	●	14Fr Male	40 cm	<b>Soft</b>	30
1064	●	14Fr Male	40 cm	Firm	30
1065	●	14Fr Male	40 cm	Firm - Curved Packaging	<b>50</b>
1082	●	16Fr Male	40 cm	Firm	30
1096	●	18Fr Male	40 cm	Firm	30
1016	●	8Fr Female	15 cm	Firm	30
1033	●	10Fr Female	15 cm	Firm	30
1047	○	12Fr Female	15 cm	Firm	30
1062	●	14Fr Female	15 cm	Firm	30

**The eyelets on all Hollister catheters are smooth and rounded.**

- Feel the difference
- Formed with ultrasonic technology
- Provide maximum comfort
- Help minimise any insertion trauma



## InstantCath: Pre-Lubricated Intermittent Catheters

- **Affordable** convenience
- **Touch free** sleeve helps protect from contamination and can be used to help direct urine flow
- **Easy and ready to use** pre-lubricated with water-soluble gel



**TOUCH FREE SLEEVE, EASY TO USE**

### InstantCath: Pre-Lubricated Intermittent Catheters

Stock Code	Funnel Colour	Size/Gender	Length	Box QTY
9670	●	10Fr Male	40 cm	25
9671	○	12Fr Male	40 cm	25
9672	●	14Fr Male	40 cm	25
9673	●	16Fr Male	40 cm	25
9674	●	10Fr Female	20 cm	25
9675	○	12Fr Female	20 cm	25
9676	●	14Fr Female	20 cm	25

## VaPro: Touch Free Hydrophilic Intermittent Catheters

It keeps it simple...

- **An innovative touch free design**
- **Designed to be evenly lubricated**
- **Helps provide trouble-free insertion and withdrawal**
- **A hydrophilic intermittent catheter**



**NEW!**  
AVAILABLE FROM AUGUST 2015

### VaPro: Touch Free Hydrophilic Intermittent Catheters

Stock Code	Funnel Colour	Size/Gender	Length	Box QTY
72122	○	12Fr Female	20 cm	30
72142	●	14Fr Female	20 cm	30
72124	○	12Fr Male	40 cm	30
72144	●	14Fr Male	40 cm	30
72164	●	16Fr Male	40 cm	30


## Advance: Pre-Lubricated Intermittent Catheters

Advance pre-lubricated intermittent catheters provide:

- **True Touch-Free Technology** – Infection risk may be reduced (environmental or anatomical)
- **Red catheter ring cap** – Keeps the Touch-Free protective tip clean and protected after the package is opened
- **Protective tip** – Helps bypass the first 15 mm of the distal urethra, where the majority of the bacteria are present. Helps reduce the risk of carrying bacteria up into the urinary system. Doubles as a guide for proper and simple insertion. See diagram below.
- **Pre-lubricated gel reservoir** – The gel amount can be user controlled
- **Smooth catheter eyelets** – Helps ensure trouble-free insertion and withdrawal



### Advance Protective Tip Technology in Action



### Advance: Pre-Lubricated Intermittent Catheters

Stock Code	Funnel Colour	Size/Gender	Length	Box QTY
92084	●	8Fr Male	40 cm	25
92104	●	10Fr Male	40 cm	25
92124	○	12Fr Male	40 cm	25
92144	●	14Fr Male	40 cm	25
92164	●	16Fr Male	40 cm	25
92184	●	18Fr Male	40 cm	25
92062	●	6Fr Female	20 cm	25
92082	●	8Fr Female	20 cm	25
92102	●	10Fr Female	20 cm	25
92122	○	12Fr Female	20 cm	25
92142	●	14Fr Female	20 cm	25