

Male Urinary Pouch



Instructions for Use

Flextend Skin Barrier



- **Extended wear time** - Median wear time of nearly 40 hours in multi-site clinical study.¹
- **Enhanced skin protection** - Resists breakdown; helps prevent skin irritation.

Latex-Free Odour-Barrier Film

- **Contains odour effectively** - Provides comfort for the patient and caregiver.
- **Prevents allergic reaction** - For latex-sensitive patients and caregivers.
- **Improves patient monitoring** - Transparent film allows easy assessment of discharge.

Application:

- Clean dry skin
- Remove pubic hair
- Measure and cut barrier
- Use ADAPT Seal if necessary
- Hold barrier in place for 30 seconds
- Check daily



The thin, tapered Flextend Skin Barrier is designed for comfort. Its superior conformability provides a secure fit for virtually any patient.

Hollister Male Urinary Pouch

Stock No.	Skin Barrier	Box QTY
9873/10	Flextend	Box of 10

For more information or samples please call Customer Engagement on:

Australia 1800 880 851
www.hollister.com.au

New Zealand 0800 678 669
www.hollister.co.nz

1. Data on file. Hollister study 1946-IC, 1997.
©2018 Hollister Incorporated. All rights reserved. Hollister and logo are trademarks of Hollister Incorporated. Printed in Australia. Oct 2019 HL072.