







CeraRings
now sealed with the
confidence of Ceramide



Adapt

Ostomy Accessories

Infused with the
science of skin health

Product	Specs/SAS Qtys	Description
Adapt CeraRing Flat Rings 	8805 Standard 8815 Slim 48mm Box of 10 SAS Monthly issue 30	<ul style="list-style-type: none"> • Infused with ceramide to help protect against dry and irritated skin • Alcohol free - No sting • Individual clam-shell packaging • Improve skin barrier fit. • Minimise leakage risk. • Protect the skin. • Improve pouch wear time. • Can be stretched, moulded or stacked to improve the fit of the skin barrier.
Adapt CeraRing Round Convex Rings 	89520 20 - 25mm 89530 30 - 35mm 89540 40 - 45mm Box of 10 SAS Monthly issue 30	<ul style="list-style-type: none"> • Alcohol free - No sting
Adapt CeraRing Oval Convex Rings 	89601 22 x 38 mm 89602 30 x 48 mm 89603 38 x 56 mm Box of 10 SAS Monthly issue 30	<ul style="list-style-type: none"> • Infused with ceramide • Alcohol free - No sting
Adapt Barrier Rings 	7806 Large 98mm Slim Box of 10 SAS Monthly issue 30	Flat, circular-shaped hydrocolloid skin barriers with a centre opening. Fill uneven skin surfaces around the stoma. Easily shaped to the dimensions of the stoma.
Adapt Barrier Extenders 	79402 20 x 25 mm Box of 30 SAS Monthly issue 60	<ul style="list-style-type: none"> • May help increase wear time • Reduces barrier edge lifting • More coverage area can enhance security • Convenient choice with up to three sections used at a time • Hydrocolloid material
Adapt Barrier Strips 	79400 Box of 10 SAS Monthly issue 30	A blend of hydrocolloids in strip format. Can be stretched and moulded to create custom shapes. Fill that narrow cavity or crease with ease.

Product	Specs/SAS Qtys	Description
Adapt Stoma Powder 	7906 28.3 gm puff bottle SAS Monthly issue 2	A blend of different hydrocolloids in a powder mixture. Absorbs excess moisture from skin around the stoma. <ul style="list-style-type: none"> • Helps protect and heal raw and weeping skin. • Helps the skin barrier adhere more effectively to most skin.
Adapt Skin Barrier Paste 	79300 57gm tube SAS Monthly issue 3 79301 14.2 gm tube Not on the SAS	A blend of hydrocolloids in a paste formulation. Fill any uneven surface. Low alcohol paste helps minimise stinging. <ul style="list-style-type: none"> • Improve skin barrier fit. • Minimise leakage risk. • Protects skin
Adapt Odour Eliminator and Lubricating Deodorant 	78500 236 ml bottle SAS Monthly issue 1 bottle 78501 8 ml sachet SAS Monthly issue 100 sachets per year (additional to bottle issue)	<ul style="list-style-type: none"> • A lubricant and deodorant in the one product. • Neutralises odour and lubricates the inside surface of the pouch to limit pancaking and ease pouch emptying
Adapt Ostomy Belt 	7300 66 - 109 cm Medium 7299 74 - 125 cm Large Each Annual issue 4	<ul style="list-style-type: none"> • An elasticised, beige, fabric belt specifically designed to attach to an ostomy pouch. • Attaches to the pouch and wraps around the body. • Helps to keep the pouch snugly in place.
Adapt Skin Gel Protective Wipes 	7917 50 sachets per box SAS Monthly issue 2 boxes	<ul style="list-style-type: none"> • Helps prevent skin stripping resulting from the frequent removal of adhesive products or skin barriers. • Helps protect skin from ostomy output.
Adapt Remover Wipes 	7760 50 sachets per box SAS Monthly issue 2 boxes	Gentle removal of adhesive products and adhesive residue. <ul style="list-style-type: none"> • Makes removal of pouching systems fast and easy. • Helps reduce the potential for skin stripping.

Product	Specs/SAS Qtys	Description
Medical Adhesive Remover Spray 	7737 76 g Spray Can Each SAS Monthly issue 2	360 degree 'Bag-in-a-can' spray technology - can be sprayed at any angle with minimal product wastage. Eases removal of skin barriers or adhesive borders (tapes), from the skin. <ul style="list-style-type: none"> • Breaks down adhesive bonds for easy removal of barrier adhesives including tapes. • Helps reduce the potential for skin stripping.
Medical Adhesive Spray 	7730 90 g Spray Can SAS Monthly issue 1	A silicone adhesive available in an aerosol can. Used to increase the adhesion of skin barriers and/or adhesive border products.* <ul style="list-style-type: none"> • Helps improve the adhesion of skin barriers. • Helps to provide a secure seal. *Requires STN authorisation
Skin Cleanser 	7210 Pump Spray Bottle 236 ml SAS Monthly issue 1	<ul style="list-style-type: none"> • Gentle foaming cleanser in a convenient pump spray. • Helps maintain healthy skin condition.
Skin Conditioning Creme 	7220 Conditioning Creme Bottle 118 ml SAS Monthly issue 2	<ul style="list-style-type: none"> • Smoothes cracked and roughened skin • Non greasy, non perfumed. • Conditions and soothes skin. • Designed for use with ostomy products. • Apply minimally.

Australia

PO Box 375
Box Hill, VIC 3128
1800 219 179

www.hollister.com.au

New Zealand

PO Box 62-027
Mt Wellington, Auckland
0800 678 669

www.hollister.co.nz

©2018 Hollister Incorporated. All rights reserved.
Hollister and logo are trademarks of Hollister Incorporated. Printed in Australia. HL707.