


A Hollister Technology

.....

AnchorFast

Oral Endotracheal Tube Fastener

*When your choice is too important
to take a chance.*



 Hollister®

NEW
AnchorFast
Slimfit™ for
Smaller Faces

Address intubated patient risks— with confidence.

Every patient dependent on a ventilator for respiration is vulnerable. Ventilator-associated risks can result in increased length of stay and increased hospital costs.¹

RISKS

of intubated patients

VENTILATOR-ASSOCIATED EVENTS such as pneumonia, due to limited access to oral care.

ORAL PRESSURE ULCERS

from the tube resting too long in one place.

UNEXPECTED TUBE OCCLUSIONS

caused by the patient.

CLINICAL BENEFITS

of the AnchorFast oral endotracheal tube fastener

EASES ACCESS TO THE ORAL CAVITY

to optimise oral care.

Findings suggest implementation of a comprehensive oral and dental care protocol can help reduce incidence of ventilator-associated pneumonia.²

HELPS PREVENT FORMATION OF LIP ULCERS

as part of a multifaceted approach.^{3,4}

HELPS PREVENT TUBE OCCLUSION

(AnchorFast Guard only) with integrated tube protection sleeve.

History, quality, and innovation.



Used at over 3,700 hospitals in the United States

In thousands of hospitals around the world, experts and professionals depend on our AnchorFast products to help improve patient care.



20 years of demonstrated quality

Hollister has a long history with endotracheal tube securement devices. Clinicians depend on its superior performance and quality, making AnchorFast products overwhelmingly the best-selling endotracheal tube securement devices on the market.




Ongoing innovation driven by customer feedback

Innovation happens when we listen to the people who use our products. What they say matters to us. It's why AnchorFast products continue to evolve to meet your needs.

AnchorFast

Oral Endotracheal Tube Fastener

The **AnchorFast** oral endotracheal tube fastener combines a number of unique features that make it the right choice.



"It's important to look for vendors who have an interest beyond just the transaction."

AnchorFast SlimFit

Oral Endotracheal Tube Fastener

The **AnchorFast SlimFit™** oral endotracheal tube fastener provides a convenient means to hold an oral endotracheal (ET) tube securely in place without the use of adhesive tape and is designed with the intent to accommodate smaller faces.

NEW
designed to
accommodate
smaller faces

NEW
upper lip
foam

Upper lip foam
designed to increase
surface area contact

Gliding tube shuttle
allows side-to-side
repositioning of the
ET tube for oral care,
while the tube remains
secure

Skin barrier pads
primary means to secure the
device to the patient's face

NEW
thumb
pad

Thumb pad with
increased surface
area for finger or
thumb placement
when opening

Grey track
that differentiates
from other
AnchorFast
products

NEW
grey
color

Adjustable neck strap device
secondary securement designed
for adjustment at an accessible
area on the patient's head

NEW
neck
strap

AnchorFast Guard

Oral Endotracheal Tube Fastener


With the addition of integrated tube protection to our trusted securement system, the **AnchorFast Guard** tube fastener can help bring a new level of confidence to the ICU.

Helps
reduce the
potential for
occlusion

Integrated tube protection sleeve





helps prevent occlusion without interfering with routine care.

Unlike separate bite blocks, it's there *before* it's needed.



AnchorFast

Oral Endotracheal Tube Fastener

	STOCK NO.	PRODUCT DESCRIPTION	QUANTITY/BOX
	9800	AnchorFast Guard Oral Endotracheal Tube Fastener	12
	9787	AnchorFast SlimFit Oral Endotracheal Tube Fastener	12
	9799	AnchorFast Oral Endotracheal Tube Fastener	12
	9788/4	AnchorFast Replacement Neck Band Straps	4

NEW
designed to
accommodate
smaller faces

Hollister Australia
P.O. Box 375,
Box Hill, Victoria, Australia 3128
www.hollister.com.au

Hollister New Zealand
58 Richard Pearse Drive,
Airport Oaks, Auckland, New Zealand
www.hollister.co.nz

For more information, please call:
Australia Customer Engagement
1800 219 179
New Zealand Customer Service
0800 678 669

¹ Restrepo M, et al. "Economic Burden of Ventilator-Associated Pneumonia Based on Total Resource Utilization." Infection Control and Hospital Epidemiology, 2010; 31(5):509-515.

² Garcia R, et al. "Reducing ventilator-associated pneumonia through advanced oral-dental care: a 48-month study." American Journal of Critical Care 18.6 (2009): 523-32.

³ Miller, Kenneth, et al. "A Utilization of AnchorFast Oral Endotracheal Tube Fastener to Reduce the Incidence of Lip Ulcers." Presented at The 55th International Respiratory Congress. 5-8 December 2009, San Antonio, Texas USA.

⁴ Veneman, Wade, et al. "An Alternative Approach to Prevent Lip Pressure Ulcers in the Intubated Patient." Presented at the Symposium on Advanced Wound Care. 17-20 April 2010, Orlando, Florida USA.