



## InView™ Silicone Male External Catheters (Male Sheaths)

Product Code	Colour	Description *	Outer Diameter	Quantity/Box
97425	Green	Special	25 mm	30
97429	Purple	Special	29 mm	30
97432	Tan	Special	32 mm	30
97436	Blue	Special	36 mm	30
97441	Orange	Special	41 mm	30
97525	Green	Standard	25 mm	30
97529	Purple	Standard	29 mm	30
97532	Tan	Standard	32 mm	30
97536	Blue	Standard	36 mm	30
97541	Orange	Standard	41 mm	30
97625	Purple	Extra	25 mm	30
97629	Purple	Extra	29 mm	30
97632	Tan	Extra	32 mm	30
97636	Blue	Extra	36 mm	30
97641	Orange	Extra	41 mm	30

### Description

- Special** Shorter Sheath Length, 50mm adhesive area
- Standard** Regular Length, 40mm adhesive area
- Extra** 70% more adhesive area than Standard



### Contact

**Australia : 1800 88 1800 880 851**  
**Email : priority@hollister.com.au**

**New Zealand : 0800 678 669**  
**Email : NZcustomerservice@libmed.co.nz**

Excerpt From AUH007 : October 2024