


A Hollister Technology

.....

AnchorFast

Oral Endotracheal Tube Fastener

*When your choice is too important
to take a chance.*



 Hollister®

AnchorFast

Oral Endotracheal Tube Fastener

The **AnchorFast** oral endotracheal tube fastener combines a number of unique features that make it the right choice.

Non-absorbent upper lip stabiliser; a new design with increased surface area is shaped to minimise pressure concentrations, keep the tube off the patient's lip and help prevent moisture buildup.

Skin barrier pads fasten the device firmly in place. The unique pad formulation is designed to be gentle on skin, easy to apply, and moisture-resistant.

Helps safeguard skin health

Gliding tube shuttle allows easy side-to-side repositioning of the ET tube for oral care, while the tube remains secure.

Helps optimise oral care

Adjustable neck strap helps provide extra security and is replaceable.

Helps stabilise ET tube

ET tube wrap with one-click security clamp locks the ET tube in place. One click releases the clamp for easy tube repositioning.

"It's important to look for vendors who have an interest beyond just the transaction. Hollister is more than a supplier, they're a partner to respiratory care practitioners in our continuous quest for advancement in patient care." — Respiratory Care Therapist

AnchorFast SlimFit

Oral Endotracheal Tube Fastener

The **AnchorFast SlimFit™** oral endotracheal tube fastener provides a convenient means to hold an oral endotracheal (ET) tube securely in place without the use of adhesive tape and is designed with the intent to accommodate smaller faces.

**Designed to
accommodate
smaller faces**

Upper lip foam

Upper lip foam
designed to increase
surface area contact

Gliding tube shuttle

allows side-to-side
repositioning of the
ET tube for oral care,
while the tube remains
secure

Skin barrier pads

primary means to secure the
device to the patient's face

Thumb pad

Thumb pad with
increased surface
area for finger or
thumb placement
when opening

Grey track
that differentiates
from other
AnchorFast
products

Grey color

Adjustable neck strap device
secondary securement designed
for adjustment at an accessible
area on the patient's head

Neck strap

"It's important to look for vendors who have an interest beyond just the transaction. Hollister is more than a supplier, they're a partner to respiratory care practitioners in our continuous quest for advancement in patient care." — Respiratory Care Therapist

AnchorFast Guard

Oral Endotracheal Tube Fastener


With the addition of integrated tube protection to our trusted securement system, the **AnchorFast Guard** tube fastener can help bring a new level of confidence to the ICU.

Helps
reduce the
potential for
occlusion

Integrated tube protection sleeve

helps prevent occlusion without
interfering with routine care.

Unlike separate bite blocks, it's there
before it's needed.



"It's important to look for vendors who have an interest beyond just the transaction. Hollister is more than a supplier, they're a partner to respiratory care practitioners in our continuous quest for advancement in patient care." — Respiratory Care Therapist

AnchorFast

Oral Endotracheal Tube Fastener

	STOCK NO.	PRODUCT DESCRIPTION	QUANTITY/BOX
	9800	AnchorFast Guard Oral Endotracheal Tube Fastener	12
	9787	AnchorFast SlimFit Oral Endotracheal Tube Fastener	12
	9799	AnchorFast Oral Endotracheal Tube Fastener	12
	9788/4	AnchorFast Replacement Neck Band Straps	4

Designed to
accommodate
smaller faces

Hollister Australia
P.O. Box 375,
Box Hill, Victoria, Australia 3128
www.hollister.com.au

Hollister New Zealand
PO Box 107097, Auckland Airport
Auckland 2150, New Zealand
www.hollister.co.nz

For more information, please call:
Australia Customer Care
1800 880 840
New Zealand Customer Service
0800 678 669

¹ Restrepo M, et al. "Economic Burden of Ventilator-Associated Pneumonia Based on Total Resource Utilization." Infection Control and Hospital Epidemiology, 2010; 31(5):509-515.

² Garcia R, et al. "Reducing ventilator-associated pneumonia through advanced oral-dental care: a 48-month study." American Journal of Critical Care 18.6 (2009): 523-32.

³ Miller, Kenneth, et al. "A Utilization of AnchorFast Oral Endotracheal Tube Fastener to Reduce the Incidence of Lip Ulcers." Presented at The 55th International Respiratory Congress. 5-8 December 2009, San Antonio, Texas USA.

⁴ Veneman, Wade, et al. "An Alternative Approach to Prevent Lip Pressure Ulcers in the Intubated Patient." Presented at the Symposium on Advanced Wound Care. 17-20 April 2010, Orlando, Florida USA.



Address intubated patient risks— with confidence.

Every patient dependent on a ventilator for respiration is vulnerable. Ventilator-associated risks can result in increased length of stay and increased hospital costs.¹

RISKS of intubated patients		CLINICAL BENEFITS of the AnchorFast oral endotracheal tube fastener	
VENTILATOR-ASSOCIATED EVENTS such as pneumonia, due to limited access to oral care.		EASES ACCESS TO THE ORAL CAVITY to optimise oral care. Findings suggest implementation of a comprehensive oral and dental care protocol can help reduce incidence of ventilator-associated pneumonia. ²	
ORAL PRESSURE ULCERS from the tube resting too long in one place.		HELPS PREVENT FORMATION OF LIP ULCERS as part of a multifaceted approach. ^{3, 4}	
UNEXPECTED TUBE OCCLUSIONS caused by the patient.		HELPS PREVENT TUBE OCCLUSION (AnchorFast Guard only) with integrated tube protection sleeve.	

History, quality, and innovation.



Used at over 3,700 hospitals in the United States

In thousands of hospitals around the world, experts and professionals depend on our AnchorFast products to help improve patient care.



20 years of demonstrated quality

Hollister has a long history with endotracheal tube securement devices. Clinicians depend on its superior performance and quality, making AnchorFast products overwhelmingly the best-selling endotracheal tube securement devices on the market.



Ongoing innovation driven by customer feedback

Innovation happens when we listen to the people who use our products. What they say matters to us. It's why AnchorFast products continue to evolve to meet your needs.